

Hardware Profiles

Hardware profiles consist of the instrument and one or more configured devices. Hardware profiles enable the software to communicate with the instrument and attached devices. You can set up these profiles using the [Hardware Configuration Editor](#), which provides an interface to configure the devices. You can set up multiple hardware profiles but only one can be active at any one time. To activate a hardware profile, the device for which the profile is configured must be connected to the workstation.

Topics in this section:

[Working with Hardware Profiles](#)

[Configuring Devices](#)

Working with Hardware Profiles

Each hardware profile must include a mass spectrometer and only peripheral devices included in the active hardware profile can be used when creating acquisition methods. Before creating an acquisition method, make sure that all devices you want to use in the method are included in the hardware profile, including a syringe pump, if your instrument comes with an integrated syringe pump. For information about troubleshooting any hardware profile issues, see [Troubleshooting hardware profiles](#).

Topics in this section:

[Creating a Hardware Profile](#)

[Adding a Peripheral Device to a Hardware Profile](#)

[Deleting a Hardware Profile](#)

[Editing an Existing Hardware Profile](#)

[About Simulation Mode](#)

Creating a Hardware Profile

1. On the Navigation bar, under **Configure**, double-click **Hardware Configuration**.
2. In the **Hardware Configuration Editor** dialog, click **New Profile**.
3. In the **Profile Name** field, type a name for the profile.
4. Click **Add Device**.
In the **Available Devices** dialog, in the **Device Type** field, **Mass Spectrometer** is the default value.
5. In the **Devices** list, select the instrument and then click **OK**.
6. In the **Devices in current profile** list, click the instrument and then click **Setup Device**.
7. (Only applicable to the 6500 series of instruments) On the **Configuration** tab, in the **Dual Mass Mode** section, select one of the following two options:
 - I **Low Mass**: To operate in limited mass range, high sensitivity operating mode, select this option. The mass range is 5–1000 Da for LIT and 50–1250 Da for RF/DC.
 - I **High Mass**: To operate in extended mass range operating mode, select this option. The mass range is 5–2000 Da for LIT and 50–2000 Da for RF/DC.
8. Select other features on the **Configuration** tab and **Communication** tab as required.
9. Click **OK** to return to the **Create New Hardware Profile** dialog.
10. Click **Add Device** and then add and configure each device that you are using with the instrument.
11. On the **Create New Hardware Profile** dialog, click **OK**.
12. To activate the hardware profile, on the **Hardware Configuration Editor**, click the hardware profile and then click **Activate Profile**.

A green check mark appears next to the profile.

Tip: You do not have to deactivate one hardware profile before activating another. Just click the hardware profile that you want to activate and then click **Activate Profile**; the other profile is deactivated automatically.

13. Click **Close**.

Note: You can view the current operating mode of the active hardware profile in the Detailed Status dialog for the instrument by double-clicking the Mass Spec icon on the lower-right section in the Analyst® software window. (only applicable to 6500 series of instruments)

Adding a Peripheral Device to a Hardware Profile

Note: Before you can add a device to an existing hardware profile, the profile must be deactivated.

1. Either [create a new hardware profile](#) or edit an existing hardware profile.
2. In the [Edit Profile](#) or the Create New Hardware Profile dialog, click **Add Device**.
The Available Devices dialog appears.
3. In the **Device Type** list, select the type of device that you want to add.
The Devices list displays a list of the available devices for the selected type.
4. In the **Devices** list, select the device that you want to add and then click **OK**.
5. In the **Devices in current profile** list, click the device and then click **Setup Device**.
6. Configure the device and then click **OK**.
7. In the Create New Hardware Profile dialog, click **OK**.
8. To activate the hardware profile, on the **Hardware Configuration Editor**, click the hardware profile and then click **Activate Profile**.
A green check mark appears next to the profile.
9. Click **Close**.

Deleting a Hardware Profile

Note: Before you can delete a hardware profile, the profile must be deactivated.

1. On the Navigation bar, double-click **Hardware Configuration**.
2. In the [Hardware Configuration Editor](#) dialog, select the hardware profile from the **Hardware Profiles** list and then click **Deactivate Profile**.
3. Click **Delete Profile** and then click **Close**.

Editing an Existing Hardware Profile

Note: Before you can add a device to an existing hardware profile, the profile must be deactivated.

1. On the Navigation bar, double-click **Hardware Configuration**.
2. In the [Hardware Configuration Editor](#) dialog, select the profile you want to edit from the list.
3. If the profile is active, then click **Deactivate Profile**.
4. Click **Edit Profile**.
5. Add, delete, or set up devices as required.
6. On the **Hardware Configuration Editor** dialog, click **OK**.

About Simulation Mode

Simulation mode is useful for creating batches and methods on a processing workstation. You can access simulation mode through the [Advanced Settings](#) dialog for each instrument or device. Each device in the active hardware configuration must be configured to run in simulation mode.

Note: The simulation options are dependent on the type of device that you are trying to simulate.

Configuring Devices

Hardware profiles contain information about the the instrument and one or more peripheral devices that you want to use with your instrument. You can configure the instrument and the devices in the Hardware Configuration Editor.

Configuring the peripheral devices requires two procedures: setting up the physical connections and configuring the software to communicate with the peripheral devices. When the software is installed, the driver required for each peripheral device is also installed. After the peripheral devices are physically connected to the computer, you can set up the appropriate configuration information.

Topics in this section:

[Configuring a Mass Spectrometer](#)

[Configuring Peripheral Devices](#)

Configuring a Mass Spectrometer

Before you can edit an existing hardware profile, you must first deactivate it.

1. Either [create a new hardware profile](#) or [edit an existing hardware profile](#).
2. If the hardware profile does not contain a mass spectrometer, click **Add Device**, and then in the **Devices** list click **Mass Spectrometer**.
3. In the **Devices in current profile** list, select the mass spectrometer and then click **Setup Device**.

The [Mass Spectrometer](#) dialog appears.

4. Select the appropriate values for the instrument.

Note: We recommend that you do not change the default settings for the GPIB communications.

If you are configuring an instrument that has an integrated Valco valve or the integrated Harvard syringe pump, select the **Use integrated injector/diverter valve** check box or the **Use integrated syringe pump** check box.

5. If you want to configure the timeouts, command retries, status poll intervals, and [simulation mode](#), on the **Communication** tab, click **Advanced** to access the [Advanced Settings](#) dialog.
6. If you want to restore the device defaults, click **Set Defaults**.
7. To save the instrument configuration, click **OK** in the Mass Spectrometer dialog.
8. To save the changes to the profile, click **OK**.

Configuring Peripheral Devices

Configuring peripherals or devices consists of two processes: setting up the physical hardware connections and configuring the software to communicate with the devices. To learn more about setting up specific peripheral devices see the *Peripheral Devices Setup Guide*.

When the software is installed, all the devices required to configure the supported devices are also installed. To activate these devices, you must configure them using the [Hardware Configuration Editor](#).