

Report Template Editor

You can print most windows, panes, or workspaces that you can see on the screen. You can also export acquisition methods, batches, quantitation Results Tables, and graph Results Tables. When exporting, you save files in formats that are appropriate for use in other programs, such as Notepad, Word, Excel, or in a Laboratory Information Management System (LIMS).

If you do not want to use the existing templates to print or export your reports, you can create new templates. You can use the Report Template Editor to customize headers, footers, and page layouts for either print or [export](#) reports. When you print a report, you send it to a printer to be printed immediately. When you export a report, you send the data to a file.

Topics in this section:

[Printing Reports](#)

[Exporting Reports](#)

[Designing Reports](#)

[Opening a Template](#)

Printing Reports

You can print reports. When you print a report, you send it to a printer to be printed immediately.

Although the process of printing varies slightly from the process of publishing, the four main steps are the same.

1. [Create a template](#) to determine the look of the report.
or
[Modify](#) an existing template.
2. [Choose](#) the report template to use for printing or exporting.
3. Select the window, pane, or workspace you want to print.
4. Print the report.

Topics in this section:

[Print Preview](#)

[Printing Without a Template](#)

[Printing with a Template](#)

[Selecting Results Table Print Options](#)

[Adjusting Print Scaling](#)

[Selecting Peak Review Sample Report Options](#)

Print Preview

You can preview, scale or copy your graphs using the **Print Preview** feature.

1. Click **File > Print Preview > Pane**.
2. Edit the Print Preview dialog.
3. Click **Print**.

Printing Without a Template

1. Click **File > Print**, and then click the item you want to print.

Printing with a Template

1. [Create a template](#) or [open](#) one that already exists.
2. Click **File > Print & Report Setup**.

The Print & Report Setup dialog appears.

3. In the **Report Template** list, click a template and then click **OK**.

Selecting Results Table Print Options

1. Open a Results Table.
2. Click **File > Print > Pane**.

From the Results Table Print Options dialog you can select Print as is or Print Grouping the Rows by Analyte.

Adjusting Print Scaling

You can adjust the size of the object you want to print.

1. Select the pane you want to scale.
2. Click **File > Print Preview > Pane**.
3. In the **Print Preview** dialog, click the **Scaling** button.
4. Adjust the scaling by a percentage or by fitting to a prescribed number of pages.

Note: The scaling button is enabled only if a pane is selected. Scaling is not supported for windows and workspaces

Selecting Peak Review Sample Report Options

1. Open peak review samples.
2. Click **File > Print > Pane**.
3. In the **Peak Review Sample Report Options** dialog, select **Print Current Sample**, or **Print all samples**.

Exporting Reports

When you export a report, you send the data to a file. You can export reports in the following formats:

- I .CSV, .DOC (Acquisition Method Editor)
- I .TXT (Batch Editor)
- I .TXT (Quantitation)
- I .txt (Graphs and Spectra)

Topics in this section:

[What Types of Reports Can I Export?](#)

[Exporting a Report](#)

What Types of Reports Can I Export?

You can export four types of reports: acquisition methods, batches, Results Tables, and graph results tables. When you export, you save files in formats that are appropriate for use in other programs, such as Notepad, Word, Excel, or in a Laboratory Information Management System (LIMS).

- I acquisition methods
- I batch details
- I graphs, spectra, or tables
- I results tables, peak reviews, or calibration curves

Exporting a Report

1. Click **File > Export**.
The **Save As** dialog appears.
2. Type the name of the file.
3. Select the file type (.CSV, .DOC, .TXT, .txt) depending on the type of report you are exporting.

4. If you are exporting a report Results Table, select either **All Columns**, or **Visible Columns** from the Export group and then click **Save**.

Designing Reports

If you are not using an existing template to print or export a report, you can create a new one or modify an existing one.

Topics in this section:

- [Creating a Report Template](#)
- [Modifying an Existing Report Template](#)
- [Adding a Design Element to Your Template](#)
- [Adding a Date](#)
- [Adding the Time](#)
- [Adding a Custom Field](#)
- [Setting Margins](#)
- [Adding a Graphic](#)
- [Adding Customized Text](#)
- [Selecting a Destination](#)
- [Using an Existing Template](#)

Creating a Report Template

If you are not using an existing template to print or export a report, you must create a new one.

1. Click **File > New > Report Template**.

The **Report Template Editor Dialog: Header/Footer Design** tab appears.

Note: Individual options are available to users depending upon security level.

2. Select a header or footer box.

The box appears outlined in red on the screen. This outline will not print on the document.

3. In the **Fields** list on the left, expand the required elements.

4. Double-click a **design element**. If the field chosen is Date, Time, or Custom Field, fill in the dialog that appears and then click **OK**.

The design element is added to the selected header or footer box.

5. Repeat steps 3 to 4 until you have completed your template. You do not need to fill in all of the boxes.

6. Click the **Page Setup** tab.

7. Select **Landscape** or **Portrait** from the Page Orientation.

8. **Set margins** to provide the appropriate space for the header and footer you have just created.

The value must be set to some value other than zero to enable a proper print out. Default margins are provided for a total of 4 lines of header or footer using the default font.

9. If required, in the **From Edge** group, type the distance from the distance from the edge of the page to the first line of the header and the last line of the footer.

10. Click **File > Save**.

Modifying an Existing Report Template

You can edit a template you have created, or you can modify an existing template to create a new one that suits your own needs. Start with a template that is similar to what you want to create.

1. On the Navigation bar, click the **Report Template Editor**.

2. Click **File > Open**. Make sure that the Report Template (*.rpt) is selected for the **Files of type** field.
3. Double-click one of the files listed in the **Files** window.
4. Modify the template as necessary. For more information, see [Adding an Element to Your Template](#).
5. [Save the file](#) with either the same or a new name, depending on whether you want to write over the original template or not.

Adding a Design Element to Your Template

You can customize your reports by adding different types of information to the headers and footers. There are three types of information that can be added to you reports.

- I [Basic Design Elements](#)
- I [Quantitation Elements](#)
- I [Acquisition Elements](#)

1. Click **File > Open**.
2. In the **Files of Type** field, select **Report Template (*.rpt)**.
3. Select a template.
4. In the **Header/Footer Design** tab of the **Report Template Editor** dialog select the area where you want the information to appear.
5. From the left window, select the elements.
6. Click **File > Save**.

Adding a Date

You can include the date in your design for report headers and footers.

1. In the [Header/Footer Design](#) tab of the **Report Template Editor** dialog, click the area in which you want the date to appear.
2. Expand **Design Elements**.
3. Double-click **Date**.
The **Date Format** dialog appears.
4. Select a format and then click **OK**. If you have selected the numerical option, select a separator too.

Note: If you choose the 2-digit date format, and the year is less than 70, then it is based on the year 2000.

Adding the Time

You can include the time of printing in your design for report headers or footers.

1. In the [Header/Footer Design](#) tab of the **Report Template Editor** dialog, click in the header or footer area where you want the time to appear.
2. Expand **Design Elements**.
3. Double-click **Time**.
The **Time Format** dialog appears.
4. Select a format and then click **OK**.

Adding a Custom Field

You can include your own text and/or graphics, such as a logo, in a report header or footer.

1. In the **Report Template Editor** [Header/Footer Design](#) tab, click in a header or footer area.
2. Expand **Design Elements**.
3. Double-click **Custom Field**.

The Custom Field Settings dialog appears.

4. Add the [text](#) or [graphics](#) you want to appear in the header or footer and then click **OK**.

Setting Margins

You can choose the margins for the page as well as the margins for the header and footer areas of your report.

1. In the Report Template Editor window, click the [Page Setup](#) tab.
2. In the Margins group, specify each margin (**Top**, **Bottom**, **Right**, **Left** header or footer).

Note: The header margin must be less than the top page margin for a header to show and the footer margin must be less than the bottom page margin for the footer to show.

Adding a Graphic

You can place a graphic into an area of the header or footer of a report. The graphic must be in bitmap (.bmp) format.

1. In the Report Template Editor window, click the [Header/Footer Design](#) tab.
2. Expand **Design Elements**.
3. Double-click **Custom Field**.
The [Custom Field Settings](#) dialog appears.
4. Click the **Graphics File** option.
5. Click the button to the right of **Graphics File** to navigate to the folder containing the graphics or [Acquisition Elements](#) file.
6. Select the file and then click **OK**.

Adding Customized Text

You can add customized text to an area in the header or footer of a report.

1. In the Report Template Editor window, click the [Header/Footer Design](#) tab.
2. Click **Design Elements**.
3. Double-click **Custom Field**.
The [Custom Field Settings](#) dialog appears.
4. Click the **Text** option.
5. Type text in the box and then click **OK**.

The boxes in the header and footer are of fixed size. Any text that is too long for the box will be truncated.

Selecting a Destination

If you do not specify a destination printer, then the Windows default printer is used.

1. On the Navigation bar, under **Configure** mode, double-click **Report Template**.
The [Report Template Editor](#) dialog appears.
2. Create a [new template](#).
or
[Open](#) an existing template.
3. Click the [Page Setup](#) tab.
4. In the **Select Destination** field, click **Printer**.

Using an Existing Template

You can use an existing template when you print.

1. Click **File > Print & Report Setup**.

The Print dialog appears.

2. In the **Report Template** group, select a template and then click **OK**.

Opening a Template

There are two ways to open an existing template. The method you choose depends upon whether you want to view the template settings, or are ready to print.

Topics in this section:

[Opening an Existing Template to View the Settings](#)

[Selecting an Existing Template when you are Ready to Print](#)

Opening an Existing Template to View the Settings

1. Click **File > Open**.
2. Make sure that the **Files of Type** field shows **Report Template *.rpt**, and then select a template from the list.
3. Click **OK**.

The Report Template Editor dialog: Header/Footer Design tab appears with the chosen template.

Selecting an Existing Template when you are Ready to Print

1. Click **File > Print & Report Setup**.
The Print & Report Setup dialog appears.
2. Click the required file from the **Report Template** group and then click **OK**.